

Restorative Options using Vertical Height as your Guide

Your patient's options increase with more vertical height




Implant Solutions®

www.solutionsforimplants.com

1-800-995-0626

Vertical Height
Total:
16mm

8mm
8mm



Restorative Options:

Full Zirconia
N.Dura™
Implant
Hybrid



Removable
Overdenture
with 2
Abutments



Denture
Teeth/Acrylic
over Titanium
Bar




Removable
Overdenture
with Bar




Vertical Height
Total:
20mm

10mm
10mm



Vertical Height
Total:
24mm

12mm
12mm



Vertical Height
Total:
28-32mm

14-16mm
14-16mm

